

Aluma™

Skin Renewal System with FACES™ technology

Aesthetic Enhancement with a Gentle Touch

Treat wrinkles, fine lines and improve skin laxity with a safe, gentle touch using the unique Functional Aspiration Controlled Electrothermal Stimulation (FACES) technology. The Aluma vacuum-assisted hand piece conforms the skin to allow energy to safely pass through the epidermis. Therapeutic Radio Frequency (RF) energy is deposited deep in the dermis creating a uniform, contained zone of heating for collagen contraction and renewal, with little to no discomfort.