

Novus® Varia™

Premier 3-color laser Photocoagulator

The Novus Varia is cleared by the FDA for the treatment of a wide range of retinal diseases. Sophisticated, yet simple, this system is designed to let the physician focus on the treatment, not the laser.

- **Versatility** – 3 independent diode-pumped solid-state lasers in one integrated package, with instant switching between wavelengths via intuitive touchscreen interface or remote control
- **Flexibility** – Wide choice of delivery devices, LaserLinks and eye safety filters; dual-fiber port enables rapid switching between delivery devices
- **Safety** – Automatic detection of delivery devices triggers automatic adjustment of laser power settings
- **Precision** – SureSpot™ optics provide sharply defined power distribution on the retina and low power density at the cornea and lens
- **Dependability** – CaseSaver™ mode enables continued treatment – even if one laser is disabled
- **Usability** – Air-cooled; plugs into 100-240 V, 50/60 Hz standard electrical outlet