

# AcuBlade™

## LARYNGEAL MICROSURGERY SYSTEM

The new AcuBlade represents a leap forward in otolaryngology. It is a revolutionary laser system that adds robotic capabilities for microsurgery. This new generation in CO2 laser accessories combines the SurgiTouch computerized laser scanning system with an operating microscope, enabling highly uniform and *automated* linear incisions and area ablations. The AcuBlade's precise capabilities enable successful treatment of vocal cord diseases while reducing the potential for injury to vocal ligaments. Clinical data shows close to a 50% reduction of operating time with dramatic simplification of procedures.